medivon.

User manual

Massage cushion

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1. Device operation

Before use, check whether the product is not damaged.
Place the device on a flat surface, such as a sofa or chair.

-  1. Press and hold the ON/OFF button for 2 seconds to turn on the device.
- 2. Press the power button to turn off/on the hot compress. The "hot compress is off" indicator will light up in blue and the "hot compress on" indicator will light up in orange. A buzzer will sound for each operation.
-  3. application of the procedure (M) to apply the type of massage cyclically, adapting it to the individual needs of the user,
- 4. Press and hold the power button for 2 seconds and the device will turn off.

2. Warnings

The device is not recommended for use by people suffering from the following diseases: People suffering from rashes or infectious skin diseases. People with diabetes, sensory disorders, malignant tumors, acute diseases, heart disease or high blood pressure. People suffering from bleeding tendencies, thrombosis, severe aneurysm or acute venous aneurysm and other symptoms of blood circulation disorders. People suffering from phlebitis, joint dysfunction, inflammation, weak bones, osteoporosis, etc. People wearing any electronic medical device such as a pacemaker. People who have pins or screws in their body or artificial joints. People with fever. Pregnant women.

3. Dangers

1. Do not use replacement parts that are not approved by the manufacturer.
2. Non-professional service personnel should not disassemble and replace parts of the leg massager without permission.
3. If the power cord is damaged, it must be replaced by appropriate professional personnel to avoid danger.
4. Do not use or store the product in a place with high humidity or high temperature.
5. After the massage is finished, unplug the power cord.
6. Do not use the product for more than 15 minutes at a time.
7. If any side effects occur, stop using the massager immediately.
8. If the temperature is too high or the product begins to smell bad or emit a burning smell, stop using it immediately and unplug the power.
9. Do not move the product by pulling the power cord to prevent the power cord from coming loose.
10. This product is not intended for persons who lack physical, sensory or mental capabilities or lack of experience and knowledge in this area (including children).
11. People who are not sensitive to heat should use the heating function carefully - preferably under the supervision of a caregiver.

4. Cleaning and care

Before cleaning, turn off the product and unplug the power cord. Wipe the surface of the product with a dry, soft cloth. Do not clean the product and accessories with corrosive detergents such as benzene, thinners. Store in a cool, dry place when not in use. Also avoid exposure to direct sunlight or humid environments, otherwise it may cause colors to fade or damage the product.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer hereby declares that the Medivon Bliss device complies with Directive 2014/35/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on restrictions on the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at the following internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Bliss