medivon.

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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Device manual

The device is ready for immediate use after purchase. Make sure the device is charged, then start the massager and select the appropriate operating mode according to your preferences. When the massager is turned off, hold down the POWER button for a few seconds to start the device.

When the device is turned on, repeated pressing of the POWER button alternates between two intelligent massage modes.

The heating button activates the heating mode. Cyclic pressing changes the heating power. The device has a 3-step temperature control. Using the suction button, we can change the suction power of the device cyclically. It is recommended to use level one the first time.

By using the intermittent suction button, we can get temporary relief and relaxation during the operation of the device. After pressing for a few seconds, the device will not suck in air.

The device is designed to be used with essential oils. Apply a thin layer on the skin and start a relaxing massage that will give you vital strength. If you notice weaker suction, check if the filter is clogged or use a new one if it needs to be replaced. The device is powered by a classic USB-C cable (5V 2A).

After the massage is completed, the device can be turned off by pressing the POWER button again.

The massager is ideal for people who need deep relaxation or improve blood circulation. With systematic use of the massage, you will reduce the symptoms of fatigue, insomnia or chronic pain of unknown origin.

The strength of the massage should be adjusted to your own preferences. The device must not be used by persons with injuries or who have undergone surgery in the last 3 months. Do not use the device in high temperatures and damp rooms. Absolutely keep away from water. Also, people with hypertension, hypotension or severe strokes cannot use the device. People with skin problems should consult a doctor before use.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after 13.08.2005.

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The manufacturer hereby declares that the Medivon Cushy device complies with Directive 2014/30 / EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Cushy