medivon.

User manual

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR

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Device manual

The massage mat is ready to work immediately after purchase. Just connect the cable to the power supply and turn on the ON / OFF button after lying comfortably on the mat. The device will automatically start working on medium intensity and program 1. Use the remote control to select your favorite massage area according to your preferences or start the heating mode. You can change the intensity of the massage at any time.

The shiatsu massage pillow can be used while massaging on the massage mat, or completely independently, all you need to do is unzip the pillow. Connect the power to the massage cushion and press the ON/OFF button located on the cushion to start it in shiatsu mode with heating. Press the button again to change the mode to massaging without heating. To turn off the pillow, hold the button for at least 3 seconds. Shiatsu massage is a very deep and effective massage, so it is not recommended to massage one area for more than 15 minutes in one session to avoid strains and injuries.

Fastening straps will allow you to lay the mat on an armchair or sofa

Warnings and Safety Instructions

Particular care should be taken by those who have recently undergone surgery or are undergoing treatment. People with hypertension, heart disease, cerebrovascular atherosclerosis and cerebral thrombosis. People with high fever, cancer and osteoporosis Women who are pregnant or menstruating. Using a pacemaker or foot implant. This product has a surface heating function and should only be used by people who are sensitive to heat. The product is not suitable for children to play with. This product can be used by children 12 years or older. When using this product, people with disabilities or people with limited behaviors should be under the supervision of their guardians. The product must not be used on children or unconscious people. Improper connection of the plug may result in electric shock. If you have any doubts about the safety of this product or the outlet is not suitable for the plug of this product, please contact customer service. If the product gets wet, it cannot be reused.

Do not use the product in a humid environment such as a bathroom or near a water source to avoid malfunction or danger to life and health. Do not use this product at a room temperature higher than 40 degrees Celsius or lower than 5 degrees Celsius, because the thickening of the lubricating oil may increase the wear of the machine transmission parts and affect the service life. Please do not use in the open or with prolonged exposure to hot sun, which may damage the varnish layer. Do not use the device under a cover where the ventilation opening may be blocked. Overheating of the device may occur. Do not carry the device by the power cord or use it as a handle. Never use the device with the air holes blocked. Never insert any objects into any openings. Do not use outdoors. Before use, long hair should be clipped so that it does not come into contact with the device. Do not use the warming function on areas where you have skin disorders, open wounds, swelling or inflammation. Improper use excludes any liability for damage. Do not use this device while driving or sleeping.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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The manufacturer hereby declares that the Medivon Nova device complies with Directive 2014/35 / EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30 / EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011 on restriction of the use of certain hazardous substances in electrical and electronic equipment. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Medivon Nova