medivon.

 **BEE  
ACTIVE**

*User Manual*

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## Acupressure mats and pillows

Forget about pain and stress - here are acupressure products that will help you fully relax and enjoy life again! Our Medivon Bee Active acupressure mats and pillows are extremely effective therapeutic tools that use the ancient art of acupressure to relieve pain, reduce muscle tension and stress.

## Tips to get started

- When using the mat or pillow for the first time, you may feel discomfort or slight pain - this is a completely natural phenomenon that will disappear after 2-3 minutes and will be replaced by blissful warmth spreading throughout the body.
- If the massage is too intense, try placing the mat on a soft surface, such as a sofa or mattress, and wearing light clothing (e.g. a cotton T-shirt) for the first few sessions.
- After the massage, you may notice spike marks and reddening of the skin on the body. This is a desirable effect of stimulating circulation and should disappear within 1 hour.

## Usage

The products are ready for immediate use after purchase. Place a mat or pillow on a flat surface and then lie down on it without clothing. Try not to make sudden movements to avoid the spikes digging into the skin too much. Lie down in a comfortable position for 15-20 minutes and try to breathe smoothly by inhaling through your nose and exhaling through your mouth.

### Warning!

The mat is safe to use, however, it requires consultation with a doctor if you have one or more of the following conditions and diseases: oncological diseases, open wounds, eczema, psoriasis, shingles, blood coagulation disorders, anticoagulant therapy, epilepsy, pregnancy, fever, thrombophlebitis, circulation problems.

Keep the product out of the reach of children. If you have moles, scars, warts, moles or an open wound on your skin, cover the area before the massage so as not to scratch the skin with sharp acupressure spikes.

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