medivon.

User manual
EYE MASSAGER

PL | EN | DE | FR | ES | PT | IT | CZ | NL
SV | SK | SL | EE | LT | LV | HU | RO | RU
NR | FI | EL | HR | DA | BG | TR | AR

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1. Construction of the device

1. Adjustment buckle
2. Elastic band
3. ON/OFF key/mode button
4. USB-C input
5. LED indicator

2. How to use

1. Before first use, please fully charge the device using the included USB cable.
2. Put the massager on your head and adjust the headband for comfortable use.
3. Press and hold the power button to turn on the product. Once activated, the button will act as a mode change function button. Each short press will change the operating mode:
 - Mode 1 Comprehensive Mode (default at power on)
 - Mode 2: Eye protection mode (vibration)
 - Mode 3: Eye protection mode (rhythmic vibration)
 - Mode 4: Sleep mode (gentle vibration)
4. Once the massage is complete, hold down the power button to turn off the device.

3. Charging and removing batteries

If you hear three consecutive beeps while using the device, it means that the built-in lithium battery has discharged. To charge it, you need to use a USB cable. While charging, the charging indicator flashes. It takes approximately 2-3 hours to fully charge. Battery life will be shortened if the battery is not used for a long time. It is recommended to charge the battery at least once a month.

The built-in lithium battery of this product contains substances that may pollute the environment. Please remove the built-in lithium battery before disposing of this product. Follow illustration 2 and the instructions:

1. Turn off or disconnect the power.
2. Remove the protective flap.
3. Remove the lithium battery.
4. The removed lithium battery should be safely disposed of at a location designated by a government or recycling facility.

4. Notes, warnings, contraindications

Please clean the product's eye shield before use. To clean, wipe the product with a slightly damp, soft cloth or towel. Do not use any corrosive substances.

People with sensitive skin or insensitive to heat should not use the product. Pregnant women, people suffering from heart disease, circulatory disorders, malignant tumors, cerebrovascular diseases and other acute diseases, as well as people under the care of a doctor, should not use the product. If you experience any discomfort or abnormal skin reaction while using the product, stop using it immediately and consult your doctor. People with eye problems should avoid using the product or use it under medical supervision. People who have undergone eye surgery and patients suffering from glaucoma, cataracts, retinal detachment and other eye diseases should not use the product. People wearing contact lenses should remove them before using the product.

Store the product in a dry and airy place, away from children. Do not expose it to direct sunlight, high temperature or moisture. This is not a toy. Children and people unable to express their needs should not use the product. Do not use the product in places with high temperature, flammable, humid or places with electromagnetic radiation. Do not disassemble the product without permission; Please make sure the product is free from impact and water to avoid failure.

5. Product specification

Name: Eye massager

Model: Horizon Travel 2

Power supply: 5V 1000mA

Lithium battery: 3.7V 700mAh

Input power: 5W

Package size: 200x91x85mm

Net weight: about 190g

Product normal operating conditions:

Ambient temperature: +5°C ~ +40°C

Relative humidity: ≤ 80%

Atmospheric pressure: 86 kPa ~ 106 kPa

Product transport and storage conditions:

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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The manufacturer declares under his sole responsibility that the Horizon Travel 2 device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as Directive 2011 of the European Parliament and of the Council /65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to the Directive of the European Parliament and 2011/65/EU as regards the list of restricted substances. The full text of the declaration of conformity is available at the internet address: www.medivon.pl/deklaracje

Importer: Armare S.A., | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Horizon Travel 2