medivon.

## Gun massager

Percussion therapy

PL | EN | DE | FR | ES | PT | IT | CZ | SK  
SL | EE | LT | LV | HU | RO | RU | NB | FI  
EL | HR | SV | DA | NL | TR | BG | AR

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## 1. Operation of the device

The device is very easy to use - one press of a button is enough to start a relaxing massage of the whole body.

After turning on the device, a short press of the power button cyclically changes the massage modes that are visible on the massager housing. Universal SMART mode is enabled by default upon first start-up. Each mode has a different application and number of strokes per minute.

The device has a built-in timer and will automatically turn off after 10 minutes of massage.

When you intend not to use the device for a long time, turn off the battery power at the bottom of the massager with the slider. This way we will avoid over-discharging the battery. Usually 2-3 hours are required to fully charge the battery. Do not use chargers other than the one provided for this purpose. While charging, the green LED will indicate a fully charged device. It is recommended not to fully discharge the device to extend the battery life.

Note: we only replace the massage heads when the device is turned off. It is enough to remove them from the device by pulling (they do not have a screw thread).

## 2. Massage

Start by choosing the massage head that suits you and your needs. Disassemble the previous one if it is attached to the massager and then attach the chosen one. Make sure that it is firmly attached and not damaged. Then start the device with the POWER button, choose your favorite massage mode.

The strength of the massage is adjustable by the user. Depending on the force with which we press the massager to the body, we will obtain different end effects. It is recommended to use the massager gently during the first massage.

### 3. Dangers

1. Do not use spare parts that are not approved by the manufacturer.
2. Unprofessional service should not disassemble and replace parts of the leg massager without permission.
3. If the power cord is damaged, it must be replaced by appropriate professional personnel to avoid a hazard.
4. Do not use or store the product in a place with high humidity or high temperature.
5. After completing the massage, disconnect the power cord.
6. Do not use the product for more than 15 minutes at a time.
7. If any side effects occur, please stop using the massager immediately.
8. If the temperature is too high or the product starts to smell bad or smells burnt, stop using it immediately and disconnect the power.
9. Do not carry the product by pulling on the power cord to prevent the power cord from loosening.
10. This product is not intended for people who have no physical, sensory or mental abilities or have no experience and knowledge in this field (including children).
11. People insensitive to heat should use the heating function carefully - preferably under the supervision of a caregiver.

### 4. Warnings

If you have any of the following conditions, consult your healthcare practitioner before use:

- Pregnant women, children or the elderly
- People who suffer from skin infections or have reddened and burned skin
- who use a pacemaker or similar implanted medical device;
- Patients with heart disease or malignant tumors, etc.
- People who are treated or who experience discomfort
- People suffering from osteoporosis
- People with high blood pressure or abnormal blood pressure

### 5. Cleaning and maintenance

Turn off the product and unplug the power cord before cleaning. Wipe the surface of the product with a dry, soft cloth. Do not clean the product and accessories with harsh detergents such as benzene, thinners, and store in a cool, dry place when not in use. Also avoid exposure to direct sunlight or humid environment, otherwise it may cause color fading or product damage.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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The manufacturer hereby declares that the Gun Pro X2 device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, and Commission Delegated Directive (EU) 2015/863 of March 31, 2015 amending Annex II to Directive 2011/65/EU. The full content of the declaration of conformity is available at the following address: [www.medivon.pl/deklaracje](http://www.medivon.pl/deklaracje)

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska  
Model: Gun Pro X2