medivon.

# User Manual

PL | EN | DE | FR | ES | PT | IT | CZ | NL  
SV | SK | SL | EE | LT | LV | HU | RO | RU  
NR | FI | EL | HR | DA | BG | TR | AR

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## Device operation

The device is suitable for immediate use after purchase. Make sure the device is charged, then put it around your neck and adjust the strap to your liking. When the massager is turned off, press and briefly hold the power button to turn it on. If the massager has not been inflated, air will be automatically pumped into the massager. The massage will start automatically once the appropriate filling is achieved. Pressing the power button cyclically changes between the two operating modes of the device. To turn off the device, hold the power button again for a few seconds. To speed up the massage, you can inflate the massager yourself through the air outlet hole. After completing the massage, to release the air, press the release in the air valve and fold the massager.

If the device is discharged, connect the device to a power source to charge the battery. While charging, the device will light up a red LED, which will turn green when the battery is fully charged. Estimated full charging time is 4-5 hours.

Do not use this device if you are pregnant, have diabetes, have a pacemaker, are recovering from surgery, suffer from migraines or neuropathy disorders. Avoid touching the product with wet hands. Frequent use may lead to abrasions. The massager must not be immersed in water. Only use replacement parts approved by the manufacturer. A non-professional service center should not disassemble or replace parts without permission. If the power cord is damaged, replacement should be carried out by professional personnel to avoid danger. Do not use or store the device in places with high humidity or temperature. Do not use the device for more than 10 minutes at a time. Stop using the device if any side effects occur. If the temperature of the device is too high, the product begins to smell bad or emit a burning smell, stop using it immediately. Do not move the product by pulling the power cord. The appliance is not intended for use by persons without appropriate physical, sensory or mental capabilities, including children. The massage area may turn red for a while but will return to normal after a while. Store in a dry and shaded place, at room temperature, out of reach of children.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment.

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The manufacturer declares under his sole responsibility that the Medivon Puffy device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as Directive 2011 of the European Parliament and of the Council /65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to the Directive of the European Parliament and 2011/65/EU as regards the list of restricted substances. The full text of the declaration of conformity is available at the internet address <https://medivon.pl/>

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska  
Model: Medivon Puffy