medivon.

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1. Operation of the device

The device is very easy to use - one press of a button is enough to start a relaxing full body massage that will allow you to rest after a hard day at work or effectively massage tense muscles after training.

After starting the device, a short press of the power button changes the massage modes cyclically, gradually increasing the frequency and power of vibrations. At first start-up, the lowest of the modes (one of five) is enabled by default. The sixth mode is the passive mode in which the massager is in standby mode.

To completely switch off the device, hold the power button for about 3 seconds.

It is recommended to fully charge the device battery before use with the 2A charger. For this purpose, use the USB-C cable included in the set. A fully charged battery will be indicated by a backlight. It is recommended not to fully discharge the device to extend the battery life.

Note: we only replace the massage heads when the device is turned off. It is enough to remove them from the device by pulling (they do not have a screw thread).

2. Massage and settings

Start by choosing the massage head that suits you and your needs. Disassemble the previous one if it is attached to the massager and then attach the chosen one. Make sure that it is firmly attached and not damaged. Then start the device with the POWER button, choose your favorite massage mode.

The strength of the massage is regulated by the user. Depending on the force with which we press the massager to the body, we will obtain different end effects. It is recommended to use the massager gently during the first massage.

To change the head angle settings, press the main button at the rotating head and then select the appropriate orientation - then release the button. The available arrangements are 0 ° 45 ° and 90 °.

The device operates at the following frequencies:

Mode 1 - 26Hz - 800rpm

Mode 2 - 33Hz - 2000rpm

Mode 3 - 40Hz - 2400rpm

Mode 4 - 46Hz - 2800rpm

Mode 5 - 53Hz - 3200rpm

3. Dangers

1. Do not use spare parts that are not approved by the manufacturer.
2. Unprofessional service should not disassemble and replace parts of the leg massager without permission.
3. If the power cord is damaged, it must be replaced by appropriate professional personnel to avoid a hazard.
4. Do not use or store the product in a place with high humidity or high temperature.
5. After completing the massage, disconnect the power cord.
6. Do not use the product for more than 15 minutes at a time.
7. If any side effects occur, please stop using the massager immediately.
8. If the temperature is too high or the product starts to smell bad or smells burnt, stop using it immediately and disconnect the power.
9. Do not carry the product by pulling on the power cord to prevent the power cord from loosening.
10. This product is not intended for people who have no physical, sensory or mental abilities or have no experience and knowledge in this field (including children).
11. People insensitive to heat should use the heating function carefully - preferably under the supervision of a caregiver.

4. Warnings

If you have any of the following conditions, consult your healthcare practitioner before use:

- Pregnant women, children or the elderly
- People who suffer from skin infections or have reddened and burned skin
- who use a pacemaker or similar implanted medical device;
- Patients with heart disease or malignant tumors, etc.
- People who are treated or who experience discomfort
- People suffering from osteoporosis
- People with high blood pressure or abnormal blood pressure

5. Cleaning and maintenance

Turn off the product and unplug the power cord before cleaning. Wipe the surface of the product with a dry, soft cloth. Do not clean the product and accessories with harsh detergents such as benzene, thinners, and store in a cool, dry place when not in use. Also avoid exposure to direct sunlight or humid environment, otherwise it may cause color fading or product damage.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment. The device was placed on the market after August 13, 2005.

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The manufacturer hereby declares that the Gun Tytan Pro device complies with Directive 2014/35/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and Directive 2011/65/EU of the European Parliament and of the Council of June 8, 2011 on restrictions on the use of certain hazardous substances in electrical and electronic equipment. The full content of the declaration of conformity is available at the following internet address: www.medivon.pl/deklaracje

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska
Model: Gun Tytan Pro